



IHE Integration Statement			
Vendor	Product Name	Version	Date
EndoSoft, LLC	EndoVault	3.0	1/29/2016

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	Set of DICOM Instances
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	CDA Wrapped Text Report
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	PDF Report
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	CDA Imaging Report with Structured Headings
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	None
ED Physician Note	Content Consumer	None
Laboratory Specimen Barcode Labeling	Label Broker	Query for labeling instructions
Laboratory Specimen Barcode Labeling	Label Broker	None
Laboratory Testing Workflow	Department System Scheduler/Order Filler	None
Laboratory Testing Workflow	Department System Scheduler/Order Filler	Report Fac-Simile For Order Group
Consistent Time	Time Client	None

Internet address for vendors IHE information : <a href="http://www.endosoft.com">http://www.endosoft.com</a>	
Links to Standards Conformance Statements for the Implementation	
HL7	<a href="http://www.endosoft.com/index.php/technical-solutions/hl7-integration.html">http://www.endosoft.com/index.php/technical-solutions/hl7-integration.html</a>
DICOM	<a href="http://www.endosoft.com/index.php/technical-solutions/dicom-integration.html">http://www.endosoft.com/index.php/technical-solutions/dicom-integration.html</a>
Link to general information on IHE : <a href="http://www.ihe.net">www.ihe.net</a>	



IHE Integration Statement			
Vendor	Product Name	Version	Date
EndoSoft, LLC	EndoVault	3.0	1/29/2016

Integration Profile	Actor	Integration Profile Option
Patient Demographics Query	Patient Demographics Consumer	None
Patient Demographics Query	Patient Demographics Consumer	Patient Demographics and Visit Query
Immunization Content	Content Consumer	None
Immunization Content	Content Consumer	View Option
Cross-Enterprise Document Sharing	Document Consumer	None
Audit Trail and Node Authentication	Secure Node	None
Sharing Lab Report	Content Consumer	None
Laboratory Specimen Barcode Labeling	Label Information Provider	None

Internet address for vendors IHE information : <a href="http://www.endosoft.com">http://www.endosoft.com</a>	
Links to Standards Conformance Statements for the Implementation	
HL7	<a href="http://www.endosoft.com/index.php/technical-solutions/hl7-integration.html">http://www.endosoft.com/index.php/technical-solutions/hl7-integration.html</a>
DICOM	<a href="http://www.endosoft.com/index.php/technical-solutions/dicom-integration.html">http://www.endosoft.com/index.php/technical-solutions/dicom-integration.html</a>
Link to general information on IHE : <a href="http://www.ihe.net">www.ihe.net</a>	