

IHE Integration Statement			
Vendor	Product Name	Version	Date
EndoSoft, LLC	EndoVault	3.1	1/27/2017

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Basic Patient Privacy Consents	Content Consumer/Content Creator	View option
Healthcare Provider Directory	Provider Information Consumer	None
Immunization Content	Content Consumer	Section Import
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Identifier Cross-Reference HL7 V3	Patient Identity Consumer	None
Scheduled Workflow	ADT Patient Registration	None
Cross-Enterprise Document Reliable Interchange	Document Source	None
Cross-Enterprise Clinical Documents Share (XDS.b)	Document Consumer/Document Source	None
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer/ Content Creator	View and Document Import
Exchange of Personal Health Record Content	Content Consumer/ Content Creator	View and Document Import
Patient Identifier Cross-Referencing for MPI	Patient Identity Consumer	None

Internet address for vendors IHE information : http://www.endosoft.com	
Links to Standards Conformance Statements for the Implementation	
HL7	https://www.endosoft.com/wp-content/uploads/2017/05/Endosoft-HL7-Interface-Specification-2016.pdf
DICOM	https://www.endosoft.com/wp-content/uploads/2017/05/EndoSoft-DICOM-Conformance-5-1-17.pdf
Link to general information on IHE : www.ihe.net	

IHE Integration Statement			
Vendor	Product Name	Version	Date
EndoSoft, LLC	EndoVault	3.1	1/27/2017

Integration Profile	Actor	Integration Profile Option
Consistent Time	Time Server/Time Client	None
Audit Trail and Node Authentication	Server Node/Secure Application	None
Document Encryption	Content Consumer/ Content Creator	None
ED Physician Note	Content Consumer	View option
Laboratory Specimen Barcode Labeling	Label Broker/ Label Information Provider	None
Laboratory Testing Workflow	Department System Scheduler/ Order Filler	None
Nursing Note	Content Consumer	None
Patient Demographics Query	Patient Demographics Consumer	None
Physician Reporting to Public Health – Cancer Registry	Content Consumer/ Content Creator	None
Sharing Lab Report	Content Consumer	None
Cross-enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Source	Text Report
Cross-Enterprise Sharing of Medical Summaries	Content Consumer/ Content Creator	None