

IHE Integration Statement			
Vendor	Product Name	Version	Date
EndoSoft, LLC	Endovault	3.1.2.0	7/2/18

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profile	Actor	Integration Profile Option
Advanced Patient Privacy Consents	Content Consumer	View Option
Advanced Patient Privacy Consents	Content Consumer	None
Audit Trail and Node Authentication	Secure Application	None
Audit Trail and Node Authentication	Secure Node	None
C-CDA 2.1 Diagnostic Imaging Report (V3)	Content Consumer	None
C-CDA 2.1 History and Physical (V3)	Content Creator	None
C-CDA 2.1 History and Physical (V3)	Content Creator	None
C-CDA 2.1 Procedure Note (V3)	Content Creator	None
Consistent Time	Time Client	None
Patient Administration Management	Patient Demographics Consumer	Merge
Patient Administration Management	Patient Demographics Consumer	None
Patient Administration Management	Patient Encounter Consumer	None
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None
Cross-Enterprise Document Media Interchange	Portable Media Creator	USB Document Storage
Cross-Enterprise Document Media Interchange	Portable Media Creator	None
Cross-Enterprise Document Media Interchange	Portable Media Importer	None

IHE Integration Statement			
Vendor	Product Name	Version	Date
EndoSoft, LLC	Endovault	3.1.2.0	7/2/18

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Media Interchange	Portable Media Importer	USB Document Storage

Internet address for vendor's IHE information: <a href="https://www.endosoft.com/wp-content/uploads/2018/09/IHE-Profile_2018.pdf">https://www.endosoft.com/wp-content/uploads/2018/09/IHE-Profile_2018.pdf</a>	
Links to Standards Conformance Statements for the Implementation	
<b>HL7</b>	<a href="#">Endosoft HL7 Interface Specification</a>
<b>DICOM</b>	<a href="#">EndoSoft DICOM Conformance Statement</a>
Link to general information on IHE : <a href="http://www.ihe.net">www.ihe.net</a>	